



Lab Resource Desk

Receiving, solvents, in-vitro reconstitution, freeze-thaw control, and batch documentation.

Handbook for qualified scientific laboratories
Review date 18 September 2026

RESEARCH USE ONLY

Not for human or veterinary use, consumption, injection, or administration.

CONTENTS**What is in this desk**

How to use these notes

01 Incoming inspection for research peptides

02 Solvent selection for lyophilized peptides

03 Reconstituting lyophilized peptides for laboratory assays

04 Aliquoting and freeze–thaw

05 How to read an Emerald Tides CoA

06 Documentation pack and product-page module

Forms index

Batch number format

07 How to read an Emerald Tides label

Intended use. These pages are general laboratory information for Research Use Only materials. They are not protocols for any catalog item, not instructions for administration, and not a substitute for an institution's SOP, chemical-hygiene plan, or IBC/IACUC requirements.

ORIENTATION**How to use these notes**

Emerald Tides AI is a research chemical supplier. It supplies selected research peptides and related compounds for scientific and laboratory environments. This desk exists so a receiving laboratory can inspect inbound batches, choose assay-compatible solvents, reconstitute lyophilized material for in-vitro work, control freeze–thaw, and keep an auditable file.

Write class-level chemistry here. Do not publish product-specific reconstitution charts, dilution calculators, or multi-puncture storage routines. Those formats leave the laboratory notebook and are read as use instructions.

On-mission topics

- Incoming inspection and three-way match of order, vial, and CoA.
- Solvent selection from charge, hydrophobicity, and assay compatibility.
- Physical reconstitution practice for in-vitro stocks.
- Single-use aliquots and freeze–thaw control.
- How to read identity versus purity on a batch-specific CoA.
- How to read an Emerald Tides batch number (example ET-003-RT).

Out of scope

- Any instruction for administration to humans or animals.
- Catalog-item “how we run this SKU” methods.
- Universal volume/mass tables that can be misread as dosing charts.

Incoming inspection for research peptides

Treat each shipment as a controlled research material. Confirm the record, the package, and the analytical file before the vial enters inventory.

Do not reconstitute a disputed vial. Solvent destroys the ability to evaluate the material as shipped.

On delivery

1. Photograph the outer carton, insulation, and any temperature indicator before opening.
2. Record carrier, tracking number, date and time received, and receiver initials.
3. Inspect for crush damage.
4. Confirm vial count, cap integrity, and that labels are legible (product identifier, batch, mass, RUO statement).

Batch number format

Emerald Tides batch numbers use the pattern ET-###-XX. Example: ET-003-RT. ET stands for Emerald Tides. The three-digit block is the production run number. The final two-letter block is the production code. Match that full string on the vial label to the CoA before the batch leaves quarantine.

Batch format: ET-###-XX (example ET-003-RT). ET = Emerald Tides. ### = production run number. XX = production code.

Reconcile three documents

Record	What to match
Purchase order / invoice	Catalog name, mass per vial, quantity, ship-to
Vial label	Batch number (ET-###-XX), product identifier, net content
Certificate of Analysis	Same batch number, same product identifier, test date and methods

Do not place a batch into general use until the CoA batch number matches the vial. File the PDF with the batch and keep the original filename.

Recommended batch file

- Invoice and packing list
- CoA, including any updated CoA if a retest is issued
- Receiving photographs and temperature notes
- Signed incoming inspection form
- Subsequent aliquot log

Quarantine versus release

Many laboratories use a two-bin system: quarantine until CoA and label match; release only after the receiving form is signed. That habit matters more than any single storage temperature.

If something is wrong

Photograph the label and carton. Keep the material at the labeled storage condition. Contact support with order number, batch, and photographs. Do not add solvent.

Solvent selection for lyophilized peptides

There is no universal solvent that both dissolves every sequence and leaves every assay untouched. Choose from peptide class, then confirm the downstream system. Always dissolve a small portion before committing the vial.

Sequence character	Typical first solvent	Laboratory notes
Basic / net positive (K, R, H rich)	Water; then 0.1% acetic acid if needed	Acid helps protonate basic residues.
Acidic / net negative (D, E rich)	Water; if poor, very dilute aqueous ammonia or bicarbonate	Avoid strong base. Neutralize before most assays.
Polar, modest hydrophobicity	Water; then dilute into assay buffer	PBS as a first solvent often reduces solubility because of salt.
Hydrophobic / aromatic-rich (L, I, V, F, W)	Minimal DMSO, DMF, or acetonitrile, then dilute	Confirm final organic percentage is tolerated by the assay.
Cys, Met, or Trp present	Degassed or low-oxygen aqueous solvents	Avoid DMSO with Cys or Met; DMSO can oxidize sulfur-containing residues.
Aggregation-prone	Urea or guanidine HCl only if the assay allows	Chaotropes are last-resort solubilizers, not routine diluents.

Order of operations

1. Calculate a concentrated stock required by the assay method. Many published in-vitro methods use stocks in the 1–5 mg/mL range; your system defines the number.
2. Add solvent down the vial wall, not as a jet onto the cake.
3. Allow the cake to wet. Gentle swirling or rolling. Avoid foaming; the air–water interface denatures and adsorbs peptide.
4. Brief bath sonication only if solids remain and the sequence is not known to be shear- or heat-sensitive.
5. If still incomplete, do not keep stacking unrelated solvents into the same vial. Volatile systems can be removed by lyophilization and a second system tried on the same solid.
6. Dilute the clear stock into assay buffer immediately before use when stability is unknown.

Assay compatibility

- Cell-based readouts. Residual DMSO, DMF, or acetonitrile is a confounder. Many groups keep final DMSO at or below 0.1% v/v and run solvent-only controls.
- Enzyme assays. pH and ionic strength of the reconstitution solvent carry into the reaction.
- LC-MS. TFA and high organic can be acceptable for analytical stocks and unsuitable for a live assay. Keep an analytical stock and an assay stock if methods conflict.
- Phosphate and other salts. Useful as working buffers; poor as the first solvent for many lyophilizates.

A usable stock is optically clear, with no film, gel, or fines. Cloudiness is incomplete dissolution or aggregation. Do not assume concentration from mass/volume if the liquid is not homogeneous. Small-molecule research reagents do not follow peptide charge rules.

Reconstituting lyophilized peptides for laboratory assays

Lyophilized peptide is the stable form. Most avoidable loss happens in the hour after solvent is added.

This is a general laboratory method. It is not a protocol for any one catalog item and not a use instruction outside controlled research.

Before opening

- Bring the closed vial toward room temperature in a desiccator or dry box so condensate does not wet the cake.
- Record batch (ET-###-XX), mass on the label, and target stock concentration in the notebook before adding liquid.
- Decide the aliquot plan now. Reconstituting an entire vial into one working bottle is how freeze–thaw damage starts.

Mass and concentration

Label mass is the starting figure. For quantitative work, confirm whether the CoA reports peptide content versus net weight. Counter-ion and residual moisture change the difference. Document the basis used (net content versus as-is) next to the batch number.

Physical handling

1. Add solvent down the glass wall.
2. Allow the cake to wet, then swirl. Do not vortex as a default.
3. Inspect against light for optical clarity.
4. Transfer to low-bind tubes if the peptide will be used at low concentration. Many sequences adsorb to untreated plastic.

Working dilutions

Prepare a concentrated stock in the solvent that actually dissolves the sequence. Prepare working dilutions in the assay buffer on the day of the experiment when stability is unknown. An internal SOP should record stock concentration and solvent, the dilution series, vehicle controls, and time-to-use after thaw. This desk does not publish universal dilution tables.

Sterility and filtration

For cell-free analytical work, cleanliness and particulate control matter more than pharmaceutical sterility. For cell culture, the laboratory's sterile-technique SOP governs. Filter only if the method requires it and the peptide is not lost on the membrane. Low-mass, hydrophobic peptides often are. If an assay requires a multi-day bench stock, the laboratory's SOP and institutional rules govern sterile technique and any antimicrobial additive. Single-use frozen aliquots are preferred.

Label immediately: product identifier, batch (ET-###-XX), stock concentration and solvent, date and time reconstituted, preparer, storage location, discard rule. Blended research materials complicate identity and purity interpretation; keep a separate batch file for each component batch.

Aliquoting and freeze–thaw

Ice is not a gentle storage medium. Freeze once per aliquot. Thaw once. Do not refreeze.

As water freezes, solutes concentrate in the unfrozen fraction, pH can drift, and ice surfaces nucleate aggregates. Each thaw repeats that insult. Apparent loss of activity in an assay is often loss of intact, soluble peptide.

Preferred hierarchy

1. Keep the unopened supplier vial lyophilized at the labeled condition — typically -20°C , desiccated, and dark; -80°C if the laboratory SOP specifies it.
2. Reconstitute only what the week's work requires, or reconstitute once and split immediately.
3. Single-use frozen aliquots at -80°C when available; otherwise in a validated -20°C freezer.
4. Working solution on ice or at bench only for the length of the assay.

How to aliquot

- Calculate aliquot volume from the assay plate map, not from convenience.
- Use small low-bind tubes or plates already accepted by the laboratory QMS.
- Flash-freeze when the SOP allows; otherwise freeze promptly in the dedicated freezer.
- Log parent batch (ET-003-RT format), stock concentration, aliquot ID, volume, freezer ID, and date.

Thawing

Thaw promptly at a controlled bench temperature, mix gently, use, and discard the remainder. Do not park a thawed aliquot in a refrigerator for the rest of the week as a default habit.

What you are protecting against

- Aggregation and loss of soluble monomer
- Oxidation of Met, Cys, and Trp during the concentrated, oxygen-available phase
- Deamidation of Asn and Gln in pH-shifted unfrozen fractions

If an assay is noisy after repeated use of one reconstituted vial, run a fresh aliquot from lyophilized material before redesigning the biology.

How to read an Emerald Tides CoA

A Certificate of Analysis is a batch-specific analytical record. It is not a release for human use and it is not FDA approval.

Question	Typical method	What it does not tell you
Is this the intended sequence / mass?	Mass spectrometry	How much related impurity is present
How much of the UV-active peak is the main species?	HPLC with a stated wavelength	Absolute peptide content if counter-ion and water are unreported
When was this batch tested?	Analysis date and laboratory name	Condition of the vial after your storage
Which production run is this vial?	Batch string ET-###-XX on label and CoA	Purity or content; the batch code is identity of the batch, not of the molecule

Lyophilized peptides often include counter-ions and residual moisture. If the CoA reports peptide content separately from net weight, use the reported basis in any mass/volume stock calculation and write that basis in the notebook.

Emerald Tides provides product-specific analytical documentation supported by independent third-party testing through JANOSHIK and BT Labs, where applicable. Independent third-party testing documents selected quality parameters. It does not constitute FDA approval, authorization, clearance, certification, or endorsement.

Documentation pack and product-page module

Print the companion forms into the batch file so a later audit can reconstruct what arrived, what was tested, and what was used in which assay.

Form	Code	Why it exists
Incoming inspection	ET-LAB-01	Package condition and three-way match of order, label, and CoA.
Batch file index	ET-LAB-02	Keeps the CoA in the batch folder instead of an inbox.
CoA reading sheet	ET-LAB-03	Forces identity, purity method, and internal spec to be written down.
Temperature excursion note	ET-LAB-04	Records what was observed and whether the laboratory accepted the batch.
Aliquot and discard log	ET-LAB-05	Shows freeze-thaw history instead of relying on memory.
Material acknowledgement	ET-LAB-06	Optional PI / lab-manager attestation of RUO institutional use.

Suggested identical block for every product page

For receiving laboratories: these materials are supplied lyophilized for Research Use Only. See the Lab Resource Desk for incoming inspection, solvent compatibility, in-vitro reconstitution practices, and freeze–thaw controls. Review the batch CoA before first use. Not for human or veterinary use.

Batch number — how to read it

Segment	Example	Meaning
ET	ET	Emerald Tides
###	003	Production run number
XX	RT	Production code
Full batch	ET-003-RT	Match this exact string on the vial and the CoA

Virginia, United States. All products are supplied strictly for laboratory research, in-vitro testing, and analytical purposes. Statements have not been evaluated by the Food and Drug Administration.

RESEARCH USE ONLY. Laboratory explainer for qualified scientific settings. Not for human or veterinary use, consumption, injection, or administration.

Answer. Read the vial face as a receiving document, not as a protocol. Confirm five things before the cap comes off: the supplier mark, the labeled mass, the production code, the Research Use Only line, and the batch string. Emerald Tides is the research chemical supplier named on the vial — not the manufacturer. The string to file is BAT: ET-003-RT — ET = Emerald Tides, 003 = production run, RT = production code.



Example vial face used in this explainer. Other masses use the same layout; only the 10MG block changes.

What each mark on the face means

On the label	What it is	How a receiving lab uses it
Flag mark	United States origin branding.	Not a regulatory approval, clearance, or FDA mark.
EMERALD TIDES™ + wave mark	Research chemical supplier named on the vial. Not a manufacturer mark.	Confirm the supplier on the purchase order, vial, and CoA.
10MG	Labeled net mass of the lyophilized solid as shipped (as-is).	Starting figure for stock math. Not peptide content and not a dose.
ET-RT (vertical code)	Production code on the face: ET + code RT.	Must agree with the last segment of the batch string.
FOR RESEARCH USE ONLY	Intended-use statement printed on the vial.	Laboratory research, in-vitro testing, and analytical work only. See ET-LAB-07.
BAT: ET-003-RT	Batch number for this specific vial.	Match this exact string to the CoA batch and ET-LAB-01 incoming inspection.

Receiving order. Photograph the closed vial so the BAT string is legible. Do not open the cap to “check the cake.” A label that cannot be read is a hold — not a reconstitution.

RESEARCH USE ONLY. Laboratory explainer for qualified scientific settings. Not for human or veterinary use, consumption, injection, or administration.

How to read BAT: ET-003-RT

The vertical block ET-RT is the production code. The unique batch lives in the BAT field. Write the full string in the batch file. ET-RT alone is not a batch number.

Segment	On this label	Meaning
ET	ET	Emerald Tides — the research chemical supplier
###	003	Production run number
XX	RT	Production code — same letters as the vertical ET-RT block
Full batch	ET-003-RT	The unique batch. This is what must match the CoA.

Three-way match at receiving

Do not release the vial from quarantine until these three strings agree. If any one disagrees, hold the batch, keep it closed, and file ET-LAB-04 if the carton is also in question.

Record	What to read	Must equal
Vial face	BAT: ET-003-RT	CoA batch number
Vial face	10MG	Label mass on the CoA / packing list
Vial face	ET-RT production code	Last segment of the CoA batch
Vial face	FOR RESEARCH USE ONLY	RUO statement on the CoA and product page
Purchase order	Catalog item and mass	What is printed on this vial

What this label does not tell you

- Identity of the sequence — that is the mass-spectrometry line on the CoA (ET-LAB-08).
- Purity % or HPLC wavelength — peak-area purity is not printed on the vial face.
- Peptide content versus counter-ion and residual water — 10MG is as-is net weight.
- Solvent, reconstitution volume, or how long a working solution may sit on the bench.
- Permission for human or veterinary use, consumption, injection, or administration.
- Who manufactured the solid. The face names the research chemical supplier, not the manufacturing site.

Short questions

Is 10MG a dose?

No. It is the labeled mass of solid in the closed vial. Stock concentration is a notebook calculation on a stated content basis.

Does the vertical ET-RT code identify the bottle by itself?

No. ET-RT is the production code. The unique batch is BAT: ET-003-RT. File the full string.

Face says ET-RT and CoA says ET-003-RT — is that a mismatch?

Not if the production code matches. The run number lives in the BAT field and on the CoA. Hold the batch if those two full strings disagree.

Does the flag or a high mass number change intended use?

No. Research Use Only is read from the vial line, the CoA, and the receiving file together. See ET-LAB-07.

Does EMERALD TIDES™ mean Emerald Tides manufactured the peptide?

No. Emerald Tides is the research chemical supplier. The mark identifies who supplied the vial, not who manufactured the solid.